



## High resolution low background CeBr<sub>3</sub> scintillators

CeBr<sub>3</sub> scintillation crystals offer an alternative to NaI(Tl) crystals for high resolution gamma spectrometry. Above an energy of 200 keV, the resolution is superior to NaI(Tl). CeBr<sub>3</sub> scintillation detectors do not suffer from the intrinsic La-138 background typical for La- halide scintillators.


|                                              |   |                           |
|----------------------------------------------|---|---------------------------|
| <b>Density</b>                               | : | 5.23 g / cc               |
| <b>Maximum emission</b>                      | : | 370 nm                    |
| <b>Decay time (typical)</b>                  | : | 18-25 ns (size dependent) |
| <b>Refractive index</b>                      | : | 2.09 (380 nm)             |
| <b>Photons yield</b>                         | : | Approx. 60.000 / MeV      |
| <b>Hygrosopic</b>                            | : | Yes                       |
| <b>Typical energy resolution @ 662 keV :</b> |   | 4 % FWHM                  |



| Energy (keV)      | Typical resolution CeBr <sub>3</sub> | Typical resolution NaI(Tl) |
|-------------------|--------------------------------------|----------------------------|
| 30 (129-I)        | 20 %                                 | 18 %                       |
| 59.5 (241-Am)     | 13 %                                 | 11 %                       |
| 122 (57-Co)       | 8 %                                  | 8.5 %                      |
| 662 (137-Cs)      | 4 %                                  | 7 %                        |
| 1332 (60-Co)      | 3 %                                  | 5.5 %                      |
| 2600 keV (Th-228) | 2.5 %                                | 4.0 %                      |

Maximum dimensions : 127 mm diameter, 152 mm high

Last revision : 03-12-2021 | Revision nr : 1 | Doc. name : CeBr<sub>3</sub> datasheet




### Intrinsic background in the Ac-227 complex

CeBr<sub>3</sub> crystals are characterized by a very small intrinsic background due to the presence of Ac-227. This results in a number of peaks between 1500 and 2200 keV. the typical count rate is 0.002 c/s/cc.

References :  
 Quarati et al. NIM A 729 (2013) pp 596-604  
 L.M. Fraile et al. NIM A 701 (2013) pp 235 - 242  
 U. Ackermann et al . NIM A 786 (2015) pp 5-11

Background spectrum  
 Inside 100 mm Lead shield



Last revision : 03-12-2021  
 Revision nr : 1  
 Doc. name : CeBr3 datasheet